IMAGE RECONSTRUCTION

UNDER VARIOUS CONDITIONS USING PARALLEL BEAMS

Emre Ataklı, Department of Electrical & Electronics Engineering, METU. e211294@metu.edu.tr

Abstract—Inverse Radon **Transform** based image reconstruction has main importance in biomedical engineering. In this study, some applications in back projection and filtered back projection (FBP) with various filters namely Ram-Lak (ramp), Cosine and Hanning filters will be shown. These techniques will be applied to a simple square image and the Shepp-Logan phantom. Best results as reconstructed image are obtained by FBP technique using Cosine filter.

Index Terms—Back-Projection, Filtering, Projection

I. INTRODUCTION

DR. Willem Roentgen has discovered X-rays in late 1895, in Wurtzburg, Germany, as a typical example of an accidental invention, seemingly, but an inspirational gift in reality. Roentgen was carrying out experiments with a Crookes tube, which is a lot common research tool at that time. When he applied large voltages to the tube to study the behavior of electrons emitted from the metal, he noticed that a piece of phosphorus substance shone. Upon this strange event, he started to try to understand what is really going on by doing a set of experiments. During these experiments he saw that the thing that causes the glow on phosphorus material can help take image of the human anatomy. After that, it had been understood that both light and X-rays are electromagnetic radiation whereas X-rays are different as they are high energy light. X-rays can penetrate through many objects because of being high energy light. But they penetrate differently through different materials according to their densities. For example, X-rays can penetrate through fat/muscle easier than bone. This is the basis in imaging the body with X-rays. After the imaging process, obtained X-ray data can be used to reconstruct medical images thanks to the fact that attenuation rate of X-rays in the body depends on tissue characteristics. [1]

In this study, some mathematical tools that are used in X-ray imaging will be implemented. Moreover, the effects of using different kind of filters, different number of beams and step sizes will be shown by comparing them.

II. THEORY

This study contains two parts as Projection and Image Reconstruction. (see Fig. 1,2, and 3)




Fig 1: Basic schematic of the term project.




Fig 2: A brief flowchart of the projection algorithm.




Fig 3: A brief flowchart of the back-projection algorithm.




Fig 4: Beam and t values for the sample image when $\theta = 0$ °

A. Mathematical background

The Radon transform and its inverse provide the mathematical foundation for reconstructing the tomographic image from projection data. [2]

B. Filters

Ramp Filter. The ramp filter is a high pass filter that does not permit low frequencies that cause blurring in the image. This type of filter is used to reduce the star artifact.

The Ramp Filter is a compensatory filter because it eliminates the star artifact resulting from simple back projection. High pass filters sharpen the edges of the image and enhance the edge information. A serious disadvantage of high pass filtering is the amplification of random noise in the image. In order to reduce the amplification of high-frequencies the ramp filter is combined with a low-pass filter.

The common method to reduce or remove random noise in a SPECT image is the application of smoothing filters. These filters are low-pass filters. In this study, two of the low pass filters were utilized: Hanning and Cosine Filters. [3]

Hanning Filter. The Hanning filter is a simple low-pass filter. [4]

Cosine Filter. This type of filter is the standard response multiplied by a cosine shape. [5]

C. Algorithm

In this study, these following steps were applied to take projection of the input image:

- Input image, step size and number of beams were specified by the user.
- Θ values were determined according to the step size and t values were determined according to the number of beams and the size of the image.
- Intersection points for all beams for all projection angles were found using the line equation.
- The points that are irrelevant to the image were removed.
- The relevant points were sorted.
- The midpoints and length of each line segment were calculated.
- The address, i.e. row and column data were found by using the size and midpoint data.
- All pixel values and corresponding distance products.
 were summed (in other words, taking integral)

Secondly, these following steps were applied to take back projection of the input image:

- Discrete Fourier Transform of the projection data was obtained.
- Projection in the frequency domain was multiplied by a desired filter in order to get rid of low frequency components in the projection data.
- Inverse Fourier Transform was obtained.
- Finally, the filtered projection data was back-projected by multiplying the distance by the projection data.
- The resultant image was shown by normalizing it.

III. RESULTS

With the help of this study, the effect of different type of filters, different number of beams and step sizes can be easily shown and compared each other. Below some of those results are included.

Python has been used for the implementation of the projection and back projection algorithms and simulation of their resultant images.

The images used for producing the projections for image reconstruction are given in the left-most side in Fig. 5 and Fig. 9 as a square image and Shepp-Logan phantom.

Results which act as quantitative evaluation measures are listed in tabular form in Table 1.

As seen from two different error measures, (in which MSE is Mean Squared Error), Cosine Filter is the best filter as compared to ramp filter and Hanning Filter.

IV. CONCLUSION

Although this study shows the effect of various filters, number of projections and step sizes on quality of the back projection image, this study can be extended with more images and filters to show their effects more clearly.

One of the limitations on the work is the language used in this work. It is Python. If C or another language similar in efficiency was used, algorithm would run fast.

Another limitation is the fact that this project and backprojection study was applied on an images, not the real human body.

Possible cause for the unexpected blurring in the unfiltered back-projection could not be found.

REFERENCES

- [1] https://sunnybrook.ca/research/content/?page=sri-groups-xray-info-3
- [2] Shahzad Ahmad Qureshi, Sikander M. Mirza, M. Arif, "Inverse Radon Transform-Based Image Reconstruction Using Various Frequency Domain Filters In Parallel Beam Transmission Tomography"
- [3] Maria Lyra and Agapi Ploussi, "Filtering in SPECT Image Reconstruction" International Journal of Biomedical Imaging, June 2011.
- [4] M. N. Salihin Yusoff and A. Zakaria, "Determination of the optimum filter for qualitative and quantitative 99mTc myocardial SPECT imaging," Iranian Journal of Radiation Research, vol. 6, no. 4, pp. 173– 181, 2009.
- $[5] \quad https://octave.sourceforge.io/image/function/rho_filter.html$

TABLE I
ERRORS AND ELAPSED TIMES UNDER
DIFFERENT CONDITIONS OF PROJECTIONS

DIFFERENT CONDITIONS OF PROJECTIONS							
	Shepp-Logan						
	Average Error	MSE	Elapsed Time				
Ramp	0.148	0.043	19.2				
Hanning	0.179	0.069	23.46				
Cosine	-0.108	0.044	18.96				
No filter	-2.953	9.500	18.18				
Square							
	Average Error	MSE	Elapsed Time				
Ramp	0.134	0.054	18.72				
Hanning	0.053	0.034	19.14				
Cosine	-0.076	0.026	18.54				
No filter	-0.952	1.061	15.3				




Fig 5: Images of original Square, reconstructed with and without filter, respectively. (Ramp filter, 180 fans and 100 detectors)




Fig 6: Images of the Square reconstructed with ramp, Hanning, Cosine filters, respectively. (180 fans and 100 detectors)




Fig 7: Images of the Square reconstructed with the step sizes 30,10, and 1 respectively. (Ramp filter, 100 detectors)




Fig 8: Images of the Square reconstructed with the number of beams 30, 100, 180 respectively. (Ramp filter, 180 fans)




Fig 9: Original Shepp-Logan, reconstructed with and without filter, respectively. (Ramp filter, 180 fans and 100 detectors)




Fig 10: Shepp-Logan image reconstructed with ramp, Cosine and Hanning filters, respectively. (180 fans and 100 detectors)




Fig 11: Images of the Shepp-Logan reconstructed with the step sizes 30,10, and 1 respectively. (Ramp filter, 100 detectors)




Fig 12: Images of the Shepp-Logan reconstructed with the number of beams 30, 100, 180 respectively. (Ramp filter, 180 fans)




Fig 13: Projections of the square image at 0, 45, 90, 135 degrees, respectively. (See Fig 4.)




Fig 14: Projections of the Shepp-Logan image at 0, 45, 90, 135 degrees, respectively. (See Fig 4.)

```
APPENDIX A:
                                                                           pro_bas = time.time()
                                                                           y_values = x_values = np.arange(-size/2, size/2+1)
import PySimpleGUI as sg
                                                                                             # determine x & y values on the
import sys, pickle, time
                                                                  image
from mplcursors import cursor
                                                                                                    np.linspace(-size/pow(2,1/2),
import scipy.io as sio
                                                                  size/pow(2,1/2),number_of_beams)
from scipy import signal
                                                                           carp = size * np.sqrt(2)
import numpy as np
                                                                           karsi uz = np.where(teta \leq 90,carp * np.cos((45-
import matplotlib.pyplot as plt
                                                                  teta) * pi/180),carp * np.cos((135-teta) * pi/180))
from numpy.fft import fft2,ifft2
                                                                           # 5. step: Find all intersection points for all beams for
from mpl_toolkits.axes_grid1 import make_axes_locatable
                                                                  all projection angles using line equation:
sg.change_look_and_feel('DefaultNoMoreNagging')
                                                                           result=[]
layout = [
                                                                           for aci in teta_degree:
                                                      # Here's
                                                                                    tan = np.tan(aci)
for the GUI window
                                                                                    \cos = \text{np.}\cos(\text{aci})
         [sg.Text('Choose where you get the projection data
                                                                                    for t degeri in t:
from:')].
                                                                                              for x degeri in x values:
  [sg.Radio('From text file
                               ', "RADIO2"), sg.Radio('From
                                                                                                      resulted_y_values = tan *
mat file ', "RADIO2"),
                                                                  x degeri + t degeri / cos #line equation
         sg.Radio('Do new projection
                                                 ', "RADIO2",
default=True)],
                                                                           result.append([aci,t_degeri,x_degeri,resulted_y_value
         [sg.Text('Enter the number of beams:')],
                                                                  s])
         [sg.InputText('100')],
                                                                           for aci in teta_degree:
         [sg.Text('Enter the step size:')],
                                                                                    \cos = \text{np.}\cos(\text{aci})
         [sg.InputText('30')],
                                                                                    \sin = \text{np.sin}(\text{aci})
         [sg.Text('kare_kosede_50ye50.mat is the default')],
                                                                                    for t_degeri in t:
         [sg.Listbox(values=['cameraman 256 256.mat','bird
                                                                                             for y_degeri in y_values:
472 472.mat',
                 'lena 256ya256.mat', 'horse 400 400.mat',
                                                                                                      if
                                                                                                                aci = 0
                                                                                                                              and
'Shepp-Logan.mat'],
                                                                  y_degeri==t_degeri:
         default_values=['kare_kosede_50ye50.mat'],
                                                                           # in case of 0 in the denominator
size=(30, 5)],
                                                                                                                for x_degeri in
         [sg.Text('Choose filter type:')],
                                                                  x_values:
         [sg.Radio('Ramp
                                 ', "RADIO3", default=True),
sg.Radio('Hanning ', "RADIO3"),
                                                                           result.append([aci,t_degeri,x_degeri,y_degeri])
         sg.Radio('Cosine ', "RADIO3"), sg.Radio('No filter
                                                                                                       elif aci != 0:
', "RADIO3")],
         [sg.Checkbox('Do
                                                                           resulted_x_values = (y_degeri * cos - t_degeri)/sin #
                                     only
                                                   projection',
default=False),sg.Text(' '*15+'Enter the projection angle:'),
                                                                  line equation
         sg.InputText(size=(5,1))],
         [sg.Submit(), sg.Cancel()]]
                                                                           result.append([aci,t_degeri,resulted_x_values,y_dege
window = sg.Window('Projection GUI', auto_size_text=True,
                                                                  ri])
default_element_size=(40, 1)).Layout(layout)
                                                                           # Remove the repeated points:
while True:
                                                                           final\_result = [list(t) for t in set(tuple(element) for
         event, values = window.Read()
                                                                  element in result)]
         if event == 'Submit':
                                                                           son = []
                  window.Close()
                                                                           # 6. Step: Remove the points which are irrelevant to
                  break
                                                                  the object:
         elif event == 'Cancel':
                                                                           # Bu işlemle irrelevant noktaları attığımız için mesela
                  window.Close()
                                                                  0 derece t=sqrt(-2) noktaları gitti
                                                                           for element in final result:
                  sys.exit()
                                                                                                                               6.5
pi = np.pi
                                                                  saniye
if values[6] == True:
                                                                                    if ((element[2]) \le (x_values[-1]) and
                                                                  (element[2]) >= (x_values[0])  and (element[3]) <= (y_values[-
         filter = 6
         filter_name = 'Ramp Filter'
                                                                   1]) and (element[3]) \geq (y_values[0])):
elif values[7] == True:
                                                                                             son.append(element)
         filter = 7
                                                                           son=sorted(son)
         filter_name = 'Hanning Filter'
                                                                                    # 7. Step: Sort the relevant points
elif values[8] == True:
                                                                  # Below, I grouped the elements of 'son' variable with respect
         filter = 8
                                                                  to their angle and t values while it had one row only before
         filter_name = 'Cosine Filter'
                                                                  this işlem
elif values[9] == True:
                                                                           temp aci t degeri = son[0][0:2]
         filter = 0
                                                                           alt_liste=[son[0]]
         filter_name = 'No Filter'
                                                                           son_son=[]
def project():
                                                                           for i in son[1:]:
```

```
if i[0:2] == temp_aci_t_degeri:
                                                                          Inputt = iter(distance_son_son)
                          alt_liste.append(i)
                                                                          son_distance_with_zeros
                           temp\_aci\_t\_degeri = i[0:2]
                                                                  [list(__import__('itertools').islice(Inputt, elem)) for elem in
                  else:
                                                                  grup]
                           son son.append(alt liste)
                                                                          # pad the projection with 0s which occur when the
                           alt liste = []
                                                                 teta values other than 45 and 90 degrees
                           alt_liste.append(i)
                                                                          son_projection_with_zeros_yeni = np.empty(0)
                           temp aci t degeri = i[0:2]
                                                                          for pro in son projection with zeros:
                                                                                   a = int( (number_of_beams - len(pro)) / 2)
         son_son.append(alt_liste)
         # 8. Find the midpoint and the length of line
                                                                                   b = number\_of\_beams - a - len(pro)
                                                                                   ekle = np.pad(pro,(a,b),'constant')
segments:
         midX=[]
                                                                                   son_projection_with_zeros_yeni
         midY=[]
                                                                  np.append(son_projection_with_zeros_yeni,ekle)
                                                                          son_projection_with_zeros
         distance_son_son=[]
         for i in son_son:
                                                                  son_projection_with_zeros_yeni.reshape(number_of_projectio
                  temp=i[0]
                                                                  ns,number of beams).tolist()
                 distance=[]
                                                                          grup_say=0
                  for j in i[1:]:
                                                                          for pro in son_distance_with_zeros:
                                                                                   if (len(pro) < number of beams):
                           temp_midX=((j[2]+temp[2])/2)
                          temp_midY = ((j[3] + temp[3])/2)
                                                                                            for
                                                                                                                              in
                           dist_temp
                                                                 range(int((number_of_beams - grup[grup_say])/2)):
                                                    pow((j[2]-
temp[2]*(j[2]-temp[2])+(j[3]-temp[3])*(j[3]-temp[3]),1/2)
                                                                                                     pro.insert(0,0)
                          midX.append(temp_midX)
                                                                                                     pro.insert(len(pro),0)
                           midY.append(temp_midY)
                                                                                   grup_say += 1
                           distance.append(dist_temp)
                                                                          with open('projection_data.txt','w') as dosya_txt:
                           temp = i
                  distance_son_son.append(distance)
                                                                          dosya_txt.write(str(number_of_projections)+'\n'+str(
                                                                 number_of_sampling_points)+'\n')
         # 9. Detect the address (row and column data) by
                                                                                   for
                                                                                                                              in
using the midpoint data.
                                                                  range(len(son_projection_with_zeros)):
         rowdata = (np.ceil(size/2 - np.floor(midY))-1)
                                                                                            dosya_txt.write(str(k+1)+'\n')
         columndata = (np.ceil(size/2 - np.floor(midX))-1)
                                                                                                                              in
         # 10. Sum all pixel value and distance products
                                                                  son_projection_with_zeros[k]:
         say = 0
                                                                                                     dosya_txt.write(str(j)+'\n')
         projection = []
                                                                          mat_array=np.array(son_projection_with_zeros)
                                                                          #list to ndarray conversion
         for i in distance_son_son:
                                                                          column array=np.array(columndata)
                  toplam=0
                                                                          row array=np.array(rowdata)
                  for j in i:
                                                                          with open('distance_list.obj','wb') as dist:
                           toplam
                                                    (i
                                                                                   pickle.dump(son_distance_with_zeros,dist)
img[int(rowdata[say])][int(columndata[say])])
                          say=say+1
                                                                          sio.savemat('projection_datas/'+values[5][0][:-
                  projection.append(toplam)
                                                                  4]+'_projection_data.mat',
         grup = []
                                                                          mdict={ 'projection':
         sa = 0
                                                                  mat_array, 'columndata':column_array, 'rowdata':row_array, 'siz
         for te in teta:
                                                                  e':size, 'original':img
                                                                                            })
                                                                          print('projection time: ',time.time() - pro_bas)
                  if (int(te) == 45 or int(te) == 135):
                          grup.append(number_of_beams)
                                                                          if values[10] == True:
                  else:
                                                                                            # If we do projection only
                          for i in range(len(t)):
                                                                          plot_projection(t,son_projection_with_zeros,number
                                   if
                                                                  _of_sampling_points,step_size)
                                             abs(t[i])
karsi_uz[sa]/2:
                                                                          return
                                            k+=1
                                                                 son_projection_with_zeros,son_distance_with_zeros,rowdata,
                                   else:
                                                                  columndata
                                            break
                                                                  plot_projection(t,projection,number_of_sampling_points,step_
                           grup.append(number_of_beams-
k*2)
                                                                 size):
                                                                          if values[11] == ":
                  sa += 1
         # açılara göre gruplu projection:
                                                                          # plot the projection of only an angle
         Inputt = iter(projection)
                                                                                   fig, axs = plt.subplots(1,4)
         son projection with zeros
                                                                                   sayyy = 0
[list(__import__('itertools').islice(Inputt, elem)) for elem in
                                                                                   for i in axs.flatten():
grup]
                                                                                            i.plot(t.round(2),projection[sayyy])
         # açılara göre gruplu distance:
                                                                                            sayyy += 1
```

```
fig.text(0.5, 0.02, 'The beam which going
                                                                           elif filter == 8:
through the image (t coordinates)',ha='center')
                                                                                   back_projection(filter_it
                 fig.text(0.1, 0.3, 'Projection value for the
                                                                  eval('signal.cosine') ))
beams',ha='center',rotation='vertical')
                                                                           elif filter == 0:
                  plt.suptitle('Projections for '+'\nNumber of
                                                                                                     # no filter
           points:
                     '+str(number of sampling points)+'\n'+'
                                                                                    back_projection()
sampling
Step size: '+str(step size))
                                                                  def back projection(getir=None):
                                                                           back pro bas = time.time()
                 plt.figure()
                 plt.imshow(img,cmap='gray')
                                                                           if getir == None:
                 plt.title('Original image')
                                                                                    getir = image_to_be_reconstructed
                                                          plot
                                                                           # Multiply the filtered projection data with the
        elif values[11] == 'all':
sinogram
                                                                  distance:
                  fig, axes = plt.subplots(1,2)
                                                                           netice = []
                                                                           for i in getir:
        axes[1].imshow(np.array(projection).T,cmap='gray')
                                                                                    o = []
                  axes[1].set ylabel('The beam which going
                                                                                    for k in i:
through the image (t coordinates)')
                                                                                            o.append(k
                  axes[1].set_xlabel('Angle')
                                                                  np.array(distance[getir.index(i)][i.index(k)]))
                  axes[1].set title('Sinogram for '+'\nNumber
                                                                                   netice.append(o)
of sampling points: '+str(number_of_sampling_points)+'\n'+'
                                                                           son_netice = []
Step size: '+str(step_size))
                                                                           for i in netice:
                  axes[0].set title('Original')
                                                                                    ara netice=[]
                  axes[0].imshow(img,cmap='gray')
                                                                                    for k in i:
        else:
                                                                                            if type(k) == np.ndarray:
                  # plot the projection of only an angle
                  cizdirilecek aci = float(values[11])
                                                                           ara netice.append(k.tolist())
                  cizdirilecek_acinin_indexi
                                                                                            else:
np.where(teta==cizdirilecek aci)[0][0]
                                                                                                     ara netice.append(k)
                  fig, axes = plt.subplots(1,2)
                                                                                    son_netice.append(ara_netice)
        axes[1].plot(t.round(2),projection[cizdirilecek acinin
                                                                           img_back = np.zeros((size,size))
_indexi],'ro')
                                                                           sav = 0
                  axes[1].set_xlabel('The beam which going
                                                                           for i in son_netice:
through the image (t coordinates)')
                                                                                    for j in i:
                  axes[1].set_ylabel('Projection value for the
                                                                                            if not j == 0:
                                                                                                     for k in j:
beams')
                  axes[1].set title('Projections
                                                           for
'+str(cizdirilecek aci)+'$^\circ$'+'\nNumber
                                             of
                                                     sampling
                                                                          img back[int(rowdata[say])][int(columndata[say])]
points: '+str(number of sampling points)+'\n'+' Step size:
                                                                  += k.real
'+str(step_size))
                                                                                                              say += 1
                  axes[0].set_title('Original')
                                                                           max_img = np.amax(img_back)
                  axes[0].imshow(img,cmap='gray')
                                                                           img_normalized = img_back / max_img
        cursor(multiple=True)
                                                                           error_img = img - img_normalized
        plt.subplots_adjust(left=0.125, bottom=0.1, right=0.9
                                                                           # if you want error_img to be included, uncomment
, top=0.85, wspace=0.4, hspace=0.2)
                                                                  related parts
        plt.show()
                                                                           img_normalized_er
                                                                                                              error_img
def ramp_filter():
                                                                  np.amax(error_img)
        t = np.linspace(0, 1, number_of_sampling_points)
                                                                           av_err = np.mean(img_normalized_er)
        return abs(abs(signal.sawtooth(2 * pi * t))-1)
                                                                           mse = np.mean(np.square(img normalized er))
def filter_it(filter_type=None,high_pass_filter=None):
                                                                           print('back
                                                                                         projection
                                                                                                      time:
                                                                                                               ',time.time()
        fft_of_projection = fft2(image_to_be_reconstructed)
                                                                  back_pro_bas)
        if high_pass_filter is None:
                                                                           print('av_err :',av_err)
                 high_pass_filter
                                                                           print('mse :',mse)
filter_type(number_of_sampling_points)
                                                                           fig,(original,back) = plt.subplots(1,2) #,error)
        filtered_fft_of_projection = fft_of_projection
                                                                           plt.subplots_adjust(left=0.125, bottom=0.1, right=0.9
high_pass_filter
                                                                  , top=0.9 , wspace=0.4, hspace=0.2)
        return ifft2(filtered_fft_of_projection).tolist()
                                                                           original.imshow(img,cmap='gray')
def filterla():
                                                                                                  im err
                                                                  error.imshow(img normalized er,cmap='gray')
        if filter == 6:
                                                                           # error related, comment/uncomment
                 back_projection(filter_it
(high pass filter=ramp filter()))
                                                                           im back
        elif filter == 7:
                                                                  back.imshow(img_normalized,cmap='gray')
                  back_projection(filter_it ( eval('np.hanning')
                                                                           divider_b = make_axes_locatable(back)
))
```

# divider_e = make_axes_locat	ke axes locatable(error)		filterla()				
#	error	re	lated,	else:			
comment/uncomment							
cax1 = divider_b.append_axes("right", size="5%", pad=0.05)			<pre># Use ready projection data (txt or mat) if values[0] == True:</pre>				
# cax2 = divider_e.append_ax	es("right",	size='	'5%",			# from txt	
pad=0.05) # error related, comme		nent		_	with	open('projection_data.txt')	as
# original.set_title('Original ima	•		1 . 1	dosya_txt:		W. 1	10
#	error	re	lated,			# data_from_txt = dosya_txt.re	ead()
comment/uncomment	:			daarra 44a.d	l:()	lines_from_txt	=
# back.set_title('Back projected #	error	ro	lated,	dosya_txt.read	nnes()	number_of_projections	_
comment/uncomment	CITOI	10	iaicu,	int(lines_from_	txt[0])	number_or_projections	=
# error.set_title('Error')				mt(mres_nom_	_tAt[O])	number_of_sampling_points	=
(, ,		#	error	int(lines_from_	txt[1])		
related, comment/uncomment				-		_to_be_reconstructed	=
fig.colorbar(im_back,cax=cax1)			image_to_be_r			
# fig.colorbar(im_err,cax=cax2)					step_si	ze = 180/number_of_projections	j
		#	error		size = 1	mat_liste[6][0][0]	
related, comment/uncomment						$ndata = mat_liste[4].tolist()[0]$	
	r_of_samp					$ta = mat_liste[5].tolist()[0]$	
"+str(number_of_sampling_points)+"		step_	_size:			mat_liste[7]	
"+str(step_size)+" "+filter_name+".png"					with of	pen('distance_list.obj','rb') as distance_list.obj','rb')	:
<pre>plt.suptitle('number_of_samplir '+str(number_of_sampling_points)+'\n'+</pre>		cton	size:	olif ve	lues[1] ==	distance = pickle.load(dist)	
'+str(step_size)+'\n'+filter_name)		siep_	_SIZE.	em va	iiues[1] —	- IIue.	
plt.show()				# fron	n mat		
prosite (()				1101		es[5] == []:	
if values[2] == True:						values[5]	=
				['kare_kosede_	50ye50.m	at']	
# If "Do new projection" is cho	sen				mat		=
if values $[5] == []$:						datas/'+values[5][0][:-	
mat	50 (1)		=	4]+'_projection			
sio.loadmat('matlar/'+'kare_kosede_50ye50.mat') # 1. step: load the default image values[5] = ['kare_kosede_50ye50.mat'] else:					ste = list(mat.values())	[2]	
					_to_be_reconstructed = mat_liste r_of_projections	[3]	
			image_to_be_reconstructed.shape[0]				
Cisc.				mage_to_6e_1		r_of_sampling_points	=
# or other image	age			number of bea		age_to_be_reconstructed.shape[1]	
mat = sio.loadmat('ma	_	es[5][0])			_to_be_reconstructed	=
img = list(mat.values())[3]				image_to_be_reconstructed.tolist()			
size = img.shape[0]					step_si	ze = 180/number_of_projections	
						mat_liste[6][0][0]	
· •	rmine the	size c	of the			$ndata = mat_liste[4].tolist()[0]$	
image	1	. C 1				$ta = mat_liste[5].tolist()[0]$	
<pre>number_of_sampling_points = int(values[3]) # 3. step:</pre>						mat_liste[7]	
int(values[3]) # 3. step: beams	gei	numbe	1 01		with of	pen('distance_list.obj','rb') as dist: distance = pickle.load(dist)	•
step_size = float(values[4])				filterl	a()	distance – pickie.ioad(dist)	
step_size = nout(varues[+])				men	u()		
# get s	tep_size			APPENDIX B	:		
$teta = np.arange(0.180, step_size)$	-						
		#		• Be sure th	at you a	re in the project code main for	older.
specify angle values according to the ste	p size					of Python are installed using	
teta_degree = teta * pi / 180			507			command prompt: "pip instal	ll −ı
number_of_projections = teta_a	idedi = teta	a.shape	[0]	requiremen			
if values[10] == True:						by writing "python projection.p	
# Do only projection						pt. Then a graphical user inte	
# Do only projection project()						r chooses the required options s, but is not interested in unre	
else:						hen after a while, resultant image	
				shown.	. 110103, 11	non arter a winie, resultant illage	oo ar C

 $image_to_be_reconstructed, distance, rowdata, column$

data = project()